

LVis 4 & tRAYcy

the basis for a new approach to gamma spectrometry.

What is LVis?

LVis is a software for **comprehensive operation and management** in a nuclear counting lab performing quantitative sample analyses!

LVis simplifies the setup and operation of spectrometers in counting laboratories and bundles extensive functionalities into a **single application**.

It uses **proven spectrum analysis techniques** and combines them with **automatically on-the-fly calculated corrections** for geometry and true coincidence summing using ORTEC's own algorithm **tRAYcy**.

ONE PLATFORM. ONE WORKFLOW. RELIABLE RESULTS.

**Full Hardware
Support**

**Automated
Measurements**

**GV Analysis
Engines**

**Automated
Calibration**

**tRAYcy MC
Modelling**

**Centralized Data
Management**

Lab Journal

**Quality
Assurance**

**Interactive Peak
Editor**

Full Traceability

GMP Compliant

**Data Import and
Export**

**Standards
Compliant**

What are LVis Key Features?

**Full Hardware
Support**

**Automated
Measurements**

**GV Analysis
Engines**

**Automated
Calibration**

**tRAYcy MC
Modelling**

**Centralized Data
Management**

Lab Journal

**Quality
Assurance**

**Interactive Peak
Editor**

Full Traceability

GMP Compliant

**Data Import and
Export**

**Standards
Compliant**

What are LVis Key Features?

What is tRAYcy?

tRAYcy is a “full Monte-Carlo” ray tracing code like MCNP/Geant4 with a clear focus on spectroscopy.

What is tRAYcy?

tRAYcy is a “full Monte-Carlo” ray tracing code like MCNP/Geant4 with a clear focus on spectroscopy.

Capable of:

- Geometry/ Attenuation correction
- Summing correction
- Mathematical calibration
- Spectrum generation
- Detector optimization

tRAYcy Detector

a digital twin in LVIS

tRAYcy is a “full Monte-Carlo” ray tracing code like MCNP/Geant4 with a clear focus on spectroscopy.

Detector details

Crystal

X-tal material: Germanium

X-tal diameter [mm]: 69.9

X-tal length [mm]: 73.4

Bulletizing radius [mm]: 8.5

Top dead layer [mm]: 1.05

Side dead layer [mm]: 1.05

X-tal hole (cavity) diameter [mm]: 8.9

X-tal hole (cavity) length [mm]: 59.7

Housing

End cap material: Aluminum

End cap diameter [mm]: 83

End cap thickness [mm]: 1

Mount cup material: Aluminum


Mount cup thickness [mm]: 0.8

Window material: Aluminum

Window-to-X-tal gap [mm]: 4

Window thickness [mm]: 1

Window diameter [mm]: 81



OK Cancel

tRAYcy hardware properties

Conversion gain: 8192

FWHM [keV] at 122 keV: 0.75 1332 keV: 1.8

Rise time [μs]: 12

Flat top time [μs]: 0.8

Gain (0.1...10): 0.6666 -> 0 - 3000 keV

OK Cancel

tRAYcy Detector

a digital twin in LVIS

Marinelli details


Container material: Polypropylene

Container diameter [mm]: 140

Container wall thickness [mm]: 1.5

Container hole diameter [mm]: 85

Container hole length [mm]: 71



OK Cancel

Reference source

Name: 1l Marinelli

Reference date: 01/07/2025

Reference time: 12:00:00

Reference volume unit: Sample

Content:

Nuclide	Activity [Bq]	Emiss./s	Energy [keV]	1-Sigma Uncert. [%]
Pb-210	16530	702.856	46.5	1.5
Am-241	4257	1529.11	59.5	1.5
Cd-109	20880	764.208	88.0	1.5
Co-57	1911	1634.1	122.1	1.5
Ce-139	2652	2118.95	165.9	1.5
Cs-137	2493	2118.8	661.7	1.5

Library: (created on 25/09/2014 09:18:20, last edit on 28/07/2025 12:32:25)

Description: Certificate No.: 1035-SE-4 0648 - 21
Type CBSS 2
Serial No.: 140921-1844008

QA library

OK Cancel

tRAYcy Detector

a digital twin in LVis


LVIS/tRAYcy live demo

tRAYcy - full Monte-Carlo ray tracing

- Scene constructed of primitives
 - cylinders/rings/spheres/..
- Ray sampling
 - homogeneous positions and directions
 - samples photon/electron/alpha cascades from level populations across decay chains
- Ray tracing (for each ray)
 - Considers photo/compton/pair interaction
 - Traces the ray until it has deposited its full energy or leaves the scene
 - Records total energy loss and energy loss in the active detector volume
 - Traces photons/electrons/positrons/alphas
 - Takes signal processing into account
- True and random summing peaks and single and double escape peaks emerge naturally




Beaker



Coax HPGe

tRAYcy Decay Data



Australian
National
University

BRICC CONVERSION COEFFICIENT
CALCULATOR




Research School of Physics
ANU College of Science

NIST

tRAYcy Decay Data

```
137BA 137CS B- DECAY 1983Be18,1996Bi23,1997WazZ07NDS 200709
137BA H TYP=ERR$AUT=J. Tuli$DAT=10-Feb-2008$COM=fixed logft values$
137BA H TYP=ERR$AUT=J. Tuli$DAT=17-Dec-2007$COM=fixed typos noted by PNPI$
137BA H TYP=FUL$AUT=E. BROWNE, J. K. TULI$CIT=NDS 108,2173 (2007)$
137BA2 H CUT=1-Oct-2006$
137BA c Measured: E(|b{+-}) (1983Be18,1978Ch22,1968Wo02,1966Hs02,1958Yo01),
137BA2c shape of the |b spectra (1983Be18,1978Ch22,1978Gr09,1969Sc23,1966Hs02),
137BA3c longitudinal polarization of |b (1975Do14), internal bremsstrahlung
137BA4c (1977We04,1975Ba20)
137BA c 1996Bi23, 1997WazZ: measured 283|g, HPGe
137BA cG 1997WazZ, 1996Bi23 looked for but did not find a transition
137BA2cG from the 661 level to the 283.5 level.
137BA cG 2007Ni04: Measured
137BA2cG |a(K)exp(127.5 keV E3 in {+134}Cs)/|a(K)exp(661.7 keV M4 in
137BAxcG {+137}Cs)=30.01 {I15}
137BA cG 2007Ya02: Measured K x-ray ratios obtained in decay vs
137BA2cG photo-ionization. Values agreed with previous measurements and theory.
137BA3cG Other: 2006Ha36
137BA cL J Adopted values
137CS P 0.0 7/2+ 30.08 Y 9 1175.63 17
137BA N 1.0 1.0 1.0
137BA PN 3
137BA L 0.0 3/2+
137BA B 1176 15.3 2 12.079 17
137BAS B EAV=416.26 8
137BA cB E from 1966Hs02, 1176 {I3} (1958Yo01), 1176 {I15} (1983Be18)
137BA cB IB$from 100-I|b(661 level).
137BA2cB Others: 5.45 {I1} (1983Be18), 5.3 {I4} (1968Wo02), 5.4 {I3} (1969Ha05),
137BA3cB 6.0 {I5} (1966Hs02), 7.6 {I8} (1958Yo01)
137BA cB Shape: |DJ=2(no) (1966Hs02,1969Sc23,1978Gr09)
137BA L 283.50 101/2+
137BA B 0.00058 8 16.61 6 2U
137BA cB IB$from I|g. 2007Se05 have recomputed logf2UT=16.45 {I12},
137BA2cB I|b=8.7|*10{+-4} {I20}, superseding their earlier value of 13.5
137BAxcB (2006SeZY)
137BAS B EAV=334.65 8
137BA G 283.5 15.8E-4 8
137BA cG E,RI$|g from 1997Waz37, %I|g average of 5.3|*10{+-4} {I14} (1996Bi23),
137BA2cG 6.1|*10{+-4} {I10} (1997Waz37)
137BA L 661.659 311/2- 2.552 M 1 M
137BA cL T$from adopted levels
137BA B 514.03 2394.7 2 9.625 2 1U
137BAS B EAV=174.32 6
137BA cB E from 1983Be18. Others: 513.89 {I26} (1978Ch22),
137BA2cB 511.63 {I84} (1968Wo02), 514 {I1} (1966Hs02), 514 {I3} (1958Yo01)
137BA cB IB$from I|g. Others: 94.6 {I3} (1968Wo02), 94.0 {I5} (1966Hs02), 92.4
137BAxcB {I8} (1958Yo01). I|b=94.6, log| {If{+1u)t}=8.5 (2006SeZY)
137BA cB
```

```
{
  "parent_element": "Cs",
  "parent_nucleons": 137,
  "parent_protons": 66,
  "parent_neutrons": 82,
  "parent": "Cs-137",
  "abundance": 0.0,
  "abundance_uncert": 0.0,
  "parent_level_energy": 0.0,
  "parent_level_energy_uncert": 0.0,
  "parent_level_j_pi": "7/2+",
  "half_life": 30.08,
  "half_life_uncert": 0.09,
  "half_life_unit": "y",
  "half_life_seconds": 948602880.0,
  "half_life_seconds_uncert": 2838239.9999999997,
  "possible_grandparents": [
    {
      "grandparent": "Xe-137",
      "decay_mode": "B-",
      "branching": 1.0,
      "branching_uncert": 0.0,
      "q_value": 4173.0,
      "q_value_uncert": 7.0
    }
  ],
  "decay_branches": [
    {
      "decay_mode": "B-",
      "branching": 1.0,
      "branching_uncert": 0.0,
      "q_value": 1175.63,
      "q_value_uncert": 0.17,
      "child_element": "Ba",
      "child_nucleons": 137,
      "child_protons": 66,
      "child_neutrons": 81,
      "child": "Ba-137",
      "x-rays": [
        {
          "gammas": [
            {
              "version": 1,
              "type": "g",
              "energy": 283.5,
              "energy_uncert": 0.1,
              "emission_probability": 0.00000059,
              "emission_probability_uncert": 8.000000000000001e-7,
              "dose": 0.00000164,
              "dose_uncert": 2.3e-7,
              "initial_level_id": 1,
              "final_level_id": 0,
              "k_shell_conversion": 0.07862500000000001,
              "k_shell_conversion_uncert": -1,
              "l_shell_conversion": 0.01234375,
              "l_shell_conversion_uncert": -1,
              "m_shell_conversion": 0.0025912749999999998,
              "m_shell_conversion_uncert": -1,
              "total_conversion": -1,
              "electron_energy": 246.06,
              "electron_energy_uncert": -1.0,
              "multipolarity": "",
              "multipole_mixing_ratio": 0.0
            },
            {
              "version": 1,
              "type": "g",
              "energy": 661.657,
              "energy_uncert": 0.003,
              "emission_probability": 0.861,
              "emission_probability_uncert": 0.002,
            }
          ]
        }
      ]
    }
  ]
}
```



tRAYcy Detector

a digital twin in LVis

LVIS/tRAYcy live demo

tRAYcy in LVis

where we are

Feature	tRAYcy status	LVis status
Efficiency transfer	implemented	implemented
TCS correction	implemented	implemented
Spectrum generation	implemented	implemented
Mathematical calibration	implemented	in testing
Well detector support	implemented	coming soon
Parameter optimization	implemented	coming soon
2.5D object support	implemented	coming soon
Full 3D scene support	implemented	coming soon
tRAYcy Analysis Engine	in testing	
Arbitrary 3D object support	coming soon	coming soon

tRAYcy in LVIs

where we are going

New decay data manager and library tool

Cs-137

Element: Caesium (Z = 55 N = 82)

Half Life (T 1/2): $3.01\text{e}+1 \pm 9.00\text{e}-2$ (years)

Decay Constant (λ): 0.0230 (1/s)

Specific Activity: 0.0230 (Bq/g)

Decay Data Updated on: January 31, 2025

Possible Parent: **Xe-137**

Decay Mode: B-

Child: **Ba-137**

Branching Ratio (%): 100.00

Q (keV): 1175.63 ± 0.17

Fluorescence Yields: K: 0.9 ± 0.004 | L: 0.11 ± 5.000


Gammas

Energy (keV) ↑↓	Intensity (%) ↑↓
661.657 ± 0.003	85.100 ± 0.200
283.5 ± 0.1	0.0005800 ± 0.0000008


X-Rays

Energy (keV) ↑↓	Intensity (%) ↑↓
32.194 ± 0.000	3.640 ± 0.100
31.817 ± 0.000	1.9900 ± 0.0500
4.47 ± 0.00	0.9100 ± 0.0400
36.378 ± 0.000	0.67200 ± 0.01800


3D Scene composer




Q & A




tRAYcy Detector Simulation



tRAYcy Detector Simulation





The detector system: MCB Properties and Detector details

Detector details

Crystal

X-tal material: Germanium

X-tal diameter [mm]: 69.9

X-tal length [mm]: 73.4

Bulletizing radius [mm]: 8.5

Top dead layer [mm]: 1.05

Side dead layer [mm]: 1.05

X-tal hole (cavity) diameter [mm]: 8.9

X-tal hole (cavity) length [mm]: 59.7

Housing

End cap material: Aluminum

End cap diameter [mm]: 83

End cap thickness [mm]: 1

Mount cup material: Aluminum

Mount cup thickness [mm]: 0.8

Window material: Aluminum

Window-to-X-tal gap [mm]: 4

Window thickness [mm]: 1

Window diameter [mm]: 81

Diagram:

OK Cancel

Add. info: CFG-SV-83
MOBIUS-ST-DET

Import Export

☐ Show tile ☒ Open on start

OK Cancel

Status	Output	Presets
Amplifier PRO	ADC	Stabilizer

the Restore

Amplifier Type

or Feedback

Optimize

Start Auto



Stop Auto

Input Polarity

(+) ○ ● (-)

File-up Rejection ☒

Close



Then define the container

Sample container

- 40 mm point source EuCo stacked
- Cerbo100
- EFS - point source on foil in ring
- EFS Czech Metrology Point Source
- EFS Point Source - 5mm dia
- Examples
- IAEA 2002 Point
- Mixed nuclide point source film
- Petri dishes
- Point Source 30 cm
- Point Source 60mm
- Point Source
- REMSPEC
- RR MOD

New sample container

New folder

Configuration

Reference Sources

Analysis libraries

Materials

Sample container

Define the container geometry

Marinelli details

Container material: Polypropylene

Container diameter [mm]: 140

Container wall thickness [mm]: 1.5

Container hole diameter [mm]: 85

Container hole length [mm]: 71

Source Material

Container Walls (all walls & bottom are assumed to be of same thickness)

Detector

OK Cancel

Sample container

Name: Marinelli 132G

Type: Marinelli

Add. info: GA-MA & Associates 1000ml
Model 132G Snap-on lid
Max height 130mm
E & Z drawing VZ-507-001

Details


OK Cancel

LVis includes more than 50 common Marinelli, vial and beaker geometries

tRAYcy Detector Simulation


- Create a realistic spectrum via Monte Carlo simulation -

Define the tRAYcy detector
in the LVIS configuration




Many detector definitions for
various ORTEC detector
models included in LVIS


Define the “hardware”
settings



Press “Start spectrum
acquisition” and then select a
reference source



OR




Create/select a routine and start a
parameterized sample count

Specify the nuclide mix

Dr. Felix Diel
+49 (0)21 59-91 36-24
Felix.Diel@Ametek.com
<https://www.ortec-online.com>

tRAYcy Detector Simulation

- Create a realistic spectrum via Monte Carlo simulation -



tRAYcy Detector Simulation

tRAYcy simulates the spectrum...

tRAYcy Comparison

Features	LabSocs ISocs	InterWinner WinnerTrack	Angle	tRAYcy/LVis
Efficiency Transfer	⊗ Not really Always need a characterization	√	√	√
Summing Correction	√	√	√	√
Spectrum Generation	⊗	⊗	⊗	√
Sourceless Calibration	√	√	⊗	Not implemented yet Release 2
Supported Isotopes	200	~200	211	> 2000
Arbitrary 3D Geometries	√	√	⊗	Not implemented yet Release 2 & 3
Detector Characterization	√	√	⊗	Not implemented yet GeDat
We own the code	⊗	⊗	⊗	√

Q & A

