

HIDEX



Hidex Sample Exchanger SEC for HpGe systems



Hidex sample exchanger for HpGe counter

- The Hidex automated high purity germanium detector is a system developed for Helsinki University Hospital Boron Neutron capture therapy laboratory.
- Installed for HUS Helsinki University Hospital back in 2018.
 - Existing HpGe detector system

Wanted focus

- Improve the workflow with sample changer






HIDEX

First we focus on your benefit,
then we create a product,
and then we do it again.

On Point. Repeat.

Hidex SEC setup



Semi-automatic sample changing unit equipped with third-party detector (and cooling unit).



Mechanical system is based on a sample rack conveyor, robotic sample moving arm and third-party detector within lead shielding.



Robotic arm picks up the sample vial and moves it on top of the detector within lead shielding.



The Sample Changer's software architecture consists of Hidex Control Platform and third-party user interface and/or measurement software.



The third-party user interface software communicates with the separate measurement software and Hidex Control Platform (conveyor and robotics control)



User defines to be measured samples by giving the user interface rack barcode ID numbers and corresponding rack position numbers.

Project based manufacturing

- To date 3 systems manufactured and taken in use
 - Neutron therapeutics Inc. in
 - Helsinki University hospital
 - Shonan Kamakura General Hospital (SKGH) in Kanagawa Prefecture, Japan.
 - Third system in US for radioisotope manufacturing
- All existing Ortec detector systems
 - GEM80P4-95
 - GEM25P4-70
- Performance information not yet available
- Initial customer feedback very positive with workflow improvement



ORTEC® / AMETEK®



Application area

Activation foil gamma spectrometry

- Neutron activation is a widely used method for the neutron fluence rate determination
- Neutron beam characterization and daily quality assurance is performed with diluted Al-Au and Al-Mn (1% weight Au/Mn) activation foils and wires Al-Mn (2.6% weight Mn)
- Activity of the foils or wires is measured with a high-purity gamma spectrometer.

Article: Accelerator-based boron neutron capture therapy facility at the Helsinki University Hospital, Liisa Porra, etc. Acta Oncologica, 61:2, 269-273, (2022)


Example performance values

- Example data from Helsinki customer
- Background
 - 3,05 cps
 - Start - Channel: 100- @ Energy: 24,63 keV
 - Stop - Channel: 8 000- @ Energy: 1 947,13 keV
 - LVis software background macro, counting time 80000s.

More information

To be investigated with current customer base

- Performance limitations
 - Lead design
 - Sample size



Thank you!

Questions?

For more information
please do not hesitate
to contact us

info@hidex.com
www.hidex.com

Tel. +358 10 843 5570