



Improving the accuracy of uncertainty calculations with Genie 4.0

Eric TISCHENBACH
Application support manager
Europe, Africa, Middle East & Asia

28 September 2023

Helsinki





MIRION
TECHNOLOGIES

Genie 4.0 Introduction

- Genie 2000's successor
Genie 4.0 provides a range of new features:

- Interactive Reports
- Improved FWHM Calibration Option
- Support for Correlations in Calibration Standards
- Automation Improvements with Python Scripting
- More Value: ISOCS, Genie-FieldPro, QA/QC
- Electronic Downloads and Licensing



- Focus on Correlation in standards efficiency calibration





Genie 2000 / Genie 4.0

Uncertainty management


- **Efficiency calibration Uncertainty**

- Certificate used
 - Rate uncertainty (combine Abundance uncertainty and **Nuclide Activity Uncertainty**)
- **Net peak area**



Genie Certificate Editor
Activity uncertainty for k=1

- Activity Uncertainty
 - **Net peak area**
 - **Efficiency**
 - Energy line Abundance
 - Nuclide Half life
 - Sample quantity
 - Random (facultative)
 - Systematic (facultative)



Genie Nuclide library Editor
• Half life uncertainty
• Abundance uncertainty



Genie 2000 /Genie 4.0

Uncertainty management

Activity

- The calculation of the uncertainty on the activity is given by the following formula

$$\sigma_C = C \cdot \sqrt{\left(\frac{\sigma_R}{100}\right)^2 + \left(\frac{\sigma_S}{S}\right)^2 + \left(\frac{\sigma_V}{V}\right)^2 + \left(\frac{\sigma_\varepsilon}{\varepsilon}\right)^2 + \left(\frac{\sigma_y}{y}\right)^2 + \left(\frac{\sigma_K}{K}\right)^2}$$

Where

$$C = \frac{S}{V \cdot \varepsilon \cdot y \cdot T_1 \cdot K_C \cdot K_W}$$

- σ_R : Random uncertainty
- σ_S : Net peak area uncertainty
- σ_V : Sample quantity uncertainty
- σ_ε : **Efficiency uncertainty**
- σ_y : peak Abundance Uncertainty
- σ_K : Nuclide Half life Uncertainty

$$C_{A_v} = \frac{\sum_{i=1}^N \frac{A_i}{\sigma_{A_i}^2}}{\sum_{i=1}^N \frac{1}{\sigma_{A_i}^2}}$$

Formula used for
Weighted mean Activity
Multi-line nuclide

- The uncertainty calculated for the nuclide average activity will be given by the formula

$$\sigma_{A_{moy}} = \sqrt{\frac{1}{\sum_{i=1}^N \frac{1}{\sigma_{A_i}^2}}}$$



Genie 2000 / Genie 4.0

Uncertainty management

ISO 11929 MDA

ISO11929 MDA formula, where measurement Uncertainty is part of MDA calculation

ISO11929 formula for decision threshold

$$C_A^* = k_{1-\alpha} \cdot \frac{W}{t_g} \sqrt{n_{b,E} + u^2_{(n_{b,E})}}$$

$$W = \frac{1}{P_E \cdot \mathcal{E}_E \cdot V \cdot f_E}$$

Where:

- C_A^* Decision threshold
- $k_{1-\alpha}$ Confidence factor
1.645 for 95%, 1.96 for 97.5%
- $n_{b,E}$ Integral of the region
- $u^2_{(n_{b,E})}$ Integral Uncertainty
- t_g Active count time

ISO11929 MDA formula

$$C_A^\# = \frac{2 \cdot C_A^* + (k^2 \cdot W) / t_g}{1 - k^2 \cdot u_{rel}^2(w)}$$

Where:

$$u_{rel}^2(w) = \left(\frac{\sigma_P}{P}\right)^2 + \left(\frac{\sigma_{\mathcal{E}\mathcal{E}}}{\mathcal{E}}\right)^2 + \left(\frac{\sigma_V}{V}\right)^2 + \left(\frac{\sigma_f}{f}\right)^2$$

- $C_A^\#$ Detection Limit (MDA)
- C_A^* Decision threshold
- $u_{rel}^2(w)$ Measurement uncertainty

- P_E Nuclide energy line abundance
- \mathcal{E}_E Energy Efficiency
- V Sample quantity
- f_E Decay correction factor



Correlations in Calibration Standards



Interactive Reports



Improved FWHM Calibration Option



Support for Correlations in Calibration Standards



Automation Improvements with Python Scripting




More Value: ISOCS, Genie-Fieldpro, QA/QC



Electronic Download and Licensing

- The uncertainty in an efficiency calibration is propagated to the uncertainty of the unknown sample activity.
- Uncertainty propagation is different for correlated values than independent values:
 - Generally additional independent measurements improve the knowledge of the detector efficiency response, driving down the final activity uncertainty.
 - Additional fully-correlated measurements *do not* improve the knowledge of the detector efficiency response, and the final activity uncertainty may be higher than with values treated independently.
- Uncertainty propagation and management is a topic of focus by regulators and the scientific community



Correlations in Calibration Standards

Interactive Reports

Improved FWHM Calibration Option

Support for Correlations in Calibration Standards

Automation Improvements with Python Scripting

More Value: ISOCS, Genie-Fieldpro, QA/QC

Electronic Download and Licensing


Definitions:

- Two values are “correlated” if there is a relationship between the measurements
- Two values are independent if there is no information gained about one measurement from the other measurement

In our application:

- In Gamma Spectroscopy and Efficiency Calibrations, the item being “measured” is the true efficiency of the detector at various energies.
- The measurement values are the calculated “efficiency triplets” from the emission lines of the calibration standard.

Isotope	Gamma-Ray Energy, keV	Half-Life, d	Activity, Bq	Flux, s ⁻¹	Uncertainty			Calibration Method**
					u _A , %	u _B , %	U, %*	
Am-241	59.5	1.580E+05	1.356E+04	4.868E+03	0.1	1.8	3.7	4π LS
Cd-109	88.0	4.614E+02	9.354E+04	3.461E+03	0.5	2.0	4.2	HPGe
Co-57	122.1	2.717E+02	2.163E+03	1.851E+03	0.4	1.7	3.5	HPGe
Ce-139	165.9	1.376E+02	3.247E+03	2.597E+03	0.4	1.8	3.7	HPGe
Cr-51	320.1	2.770E+01	9.352E+04	9.268E+03	0.1	1.8	3.7	IC
Sn-113	391.7	1.151E+02	5.726E+03	3.720E+03	0.4	2.0	4.1	HPGe
Sr-85	514.0	6.485E+01	7.275E+03	7.166E+03	0.1	1.8	3.7	IC
Cs-137	661.7	1.099E+04	2.735E+03	2.328E+03	0.7	1.9	4.1	HPGe
Mn-54	834.8	3.121E+02	9.815E+03	9.813E+03	0.1	1.7	3.5	IC
Y-88	898.0	1.066E+02	9.258E+03	8.675E+03	0.7	1.7	3.7	HPGe
Y-88	1836.1			9.184E+03	0.7	1.7	3.7	
Zn-65	1115.5	2.439E+02	1.949E+04	9.780E+03	0.1	1.8	3.7	IC
Co-60	1173.2	1.925E+03	4.325E+03	4.319E+03	0.7	1.8	3.9	HPGe
Co-60	1332.5			4.324E+03	0.7	1.8	3.9	





Genie 4.0

Correlation improvement

Introduction

- For high precision gamma spectrometry, the uncertainty from efficiency calibration is often dominating
- Quantities in gamma spectrometry can be determined using multiple measurements
 - Efficiency as a function of Energy, weighted mean activity of radionuclides
- When these measurements are independent the uncertainty is reduced when more measurements are added
- This leads to uncertainties from the efficiency fit to be lower than the uncertainty of the individual points unless correlations between data points are accounted for





Genie 4.0 Correlation

- Correlations are relations in data that varies together in a way not expected from chance alone
- The emission rate for lines from a multiline radionuclide
- Changes the least squares fit



$$\chi^2 = \sum_{i=1}^n \frac{(y_i - f(x))^2}{\sigma_{y_i}^2} \rightarrow \chi^2 = (y - f(x))^T V^{-1} (y - f(x))$$

- This has consequences for efficiency calibration, weighted mean activity and interference correction of activities*
- But also for MDA ISO11929 calculation



Example 1 - Single radionuclide calibration



- Eu-152 calibration
 - Point source at 30 cm
 - 20% rel eff p-type detector
 - Multiple emission lines with correlation from the radionuclide activity
- Activity has a relative uncertainty of 3.1% at 99% CL



Example 2 - Mixed gamma source

- Multi-radionuclide calibration source
- The level of correlation depends on the assumptions
 - All points independent
 - Points from same radionuclide correlated
 - Points using the same calibration method correlated
 - All points correlated
- For the most conservative assumption the relative uncertainty will be slightly above the lowest relative uncertainty from the certificate.

Isotope	Gamma-Ray Energy, keV	Half-Life, d	Activity, Bq	Flux, s ⁻¹	Uncertainty			Calibration Method**
					u _A , %	u _B , %	U _i , %*	
Am-241	59.5	1.580E+05	1.358E+04	4.875E+03	0.1	1.8	3.7	4π LS
Cd-109	88.0	4.614E+02	9.367E+04	3.466E+03	0.5	2.0	4.2	HPGe
Co-57	122.1	2.717E+02	2.165E+03	1.854E+03	0.4	1.7	3.5	HPGe
Ce-139	165.9	1.376E+02	3.251E+03	2.601E+03	0.4	1.8	3.7	HPGe
Cr-51	320.1	2.770E+01	9.365E+04	9.281E+03	0.1	1.8	3.7	IC
Sn-113	391.7	1.151E+02	5.733E+03	3.725E+03	0.4	2.0	4.1	HPGe
Sr-85	514.0	6.485E+01	7.285E+03	7.175E+03	0.1	1.8	3.7	IC
Cs-137	661.7	1.099E+04	2.739E+03	2.331E+03	0.7	1.9	4.1	HPGe
Mn-54	834.8	3.121E+02	9.828E+03	9.826E+03	0.1	1.7	3.5	IC
Y-88	898.0	1.066E+02	9.270E+03	8.686E+03	0.7	1.7	3.7	HPGe
Y-88	1836.1	—	—	9.196E+03	0.7	1.7	3.7	—
Zn-65	1115.5	2.439E+02	1.951E+04	9.763E+03	0.1	1.8	3.7	IC
Co-60	1173.2	1.925E+03	4.331E+03	4.324E+03	0.7	1.8	3.9	HPGe
Co-60	1332.5	—	—	4.330E+03	0.7	1.8	3.9	—



Weighted mean activity

- Independent weighted mean

$$A = \frac{\sum \frac{A_i}{\sigma_{A_i}^2}}{\sum \frac{1}{\sigma_{A_i}^2}}, \quad \sigma_A = \sqrt{\frac{1}{\sum \frac{1}{\sigma_{A_i}^2}}}$$

- New Correlated weighted mean algo

$$A = \frac{\sum_{i,j} V_{ij}^{-1} A_j}{\sum_{i,j} V_{ij}^{-1}}, \quad \sigma_A = \sqrt{\frac{1}{\sum_{i,j} V_{ij}^{-1}}}$$

From covariance matrix $V = \begin{pmatrix} \frac{\sigma_{\varepsilon_1}^2}{\varepsilon_1^2} & \dots & V_{1n} \\ \dots & \dots & \dots \\ V_{n1} & \dots & \frac{\sigma_{\varepsilon_n}^2}{\varepsilon_n^2} \end{pmatrix}$

Energy (keV)	Eff:Corr WtM:Corr	Eff:Corr WtM:Indep	Eff:Indep WtM:Corr	Eff:Indep WtM:Indep
121.8	1.39 %	1.39 %	1.39 %	1.39 %
244.7	1.52 %	1.52 %	1.42 %	1.42 %
344.3	1.41 %	1.41 %	1.12 %	1.12 %
444.0	2.47 %	2.47 %	2.43 %	2.43 %
778.9	1.62 %	1.62 %	1.35 %	1.35 %
867.4	3.33 %	3.33 %	3.18 %	3.18 %
964.1	1.48 %	1.48 %	1.08 %	1.08 %
1085.9	1.62 %	1.62 %	1.27 %	1.27 %
1112.1	1.52 %	1.52 %	1.15 %	1.15 %
1408.0	1.63 %	1.63 %	1.61 %	1.61 %
Wt Mean	1.25 %	0.52 %	0.55 %	0.44 %

Relative uncertainty of Eu-152 source 1.20 %



Correlations are supported in Genie 4.0

- **Primary User Focus:** Updated Certificate File Editor (see right).

This is where the user defines how the calibration source standard emission lines are correlated.

- **Transparent for the user**

Genie algorithms updated for efficiency calibration fits and nuclide identification with interference correction engine (now called w/ Correlations) “NID

Nuclide	Activity (Bq)	%Uncert	Energy (keV)	Intensity (per 100 decays)	Uncertainty (per 100)	Rate (\$^{-1}\$ unit \$^{-1}\$)	%Uncert	Half-life	Uncertainty	Current Rate (\$^{-1}\$ unit \$^{-1}\$)
+ CD-109	14247.100	3.023	88.034	3.6100	0.1000	514.3203	4.1000	462.90 D	2.0000	1.95e-08
CO-57	1068.310	1.700	122.061	85.9000	1.2000	917.6783	2.2000	271.80 D	0.0500	1.65e-15
CE-139	690.209	4.000	165.853	79.8860	0.0150	551.3803	4.0000	137.64 D	0.0000	4.98e-33
HG-203	2134.180	2.495	279.197	81.4600	0.1300	1738.5029	2.5000	46.61 D	0.0010	5.74e-101
SN-113	2580.310	3.643	391.688	64.0000	2.0000	1651.3986	4.8000	115.09 D	0.0400	2.03e-39
SR-85	2866.250	2.159	514.007	96.0000	2.0000	2751.6001	3.0000	64.84 D	0.0200	1.12e-71
CS-137	1671.280	2.999	661.660	85.2100	0.0700	1424.0977	3.0000	30.00 Y	0.2000	516.8092
Y-88	10918.900	3.987	898.042	93.7000	0.3000	10231	3.9985	106.65 D	0.0400	6.07e-42



Conclusion

- Taking correlations into account for efficiency calibration and activity calculations in gamma spectrometry can increase the uncertainty of the radionuclide activities when the efficiency data points are correlated.
- Uncertainties in the independent variable, intensities, are necessary when performing interference correction when the relative uncertainty in the decay data is comparable or large compared to other uncertainties.
- Nuclide Activity Uncertainty Results may increase by 1-3%

Included in Genie 4.0



Genie 4.0 webinar presentation available on our website

[Latest in Gamma Spectroscopy: New Genie™ 4.0 Software \(on24.com\)](#)

Thank you

Special thank to our scientist team
Henrik Persson, Troy Anderson, John M Kirkpatrick, Kara Phillips

