



**Measuring mass attenuation  
coefficients for materials with  
unknown composition  
by performing transmission  
measurements with a HPGe  
detector for X-rays  
and low-energy gamma rays**

Leen Verheyen

leen.verheyen@sckcen.be



Laboratory for gammaspectrometry

Sample types - Gamma lms – efftran – Genie2K

Materials and vials

Calibration curves with transmission measurements

Implementation in routine lab

Conclusion

# Laboratory for gamma spectrometry

20 HPGe det >> 5 Low Energy, 1 Well  
typen, LEGE and 12 coaxial

- Home made gamma lins
- Automatically N2 filling system
- Weighing balance with online registration
- Digital signal analyzer LYNX (Mirion)



# Samples

## Different geometries

- +/- 3500 samples per year
- Weekly QA and monthly Background
- Measurement in different geometries
- Volumes from 4 - 2400 ml
- Filters
- Cartridges
- ...





# Gamma lims

**ORDER ID** **20386** Confirmation Date   
**C\_LIMS\_ID:**

---

**Reception date**  **Order\_status:** ■ **Firm / CC or WBS**   
**Reference**  **Contact name**   
**Order Comment:**  **CC or WBS**

**O/0/2230020386**

 **Print Order Barcode**   
  **Print LAB Sheet**

---

**Sample ID** **46304** 1 of 1 **Sample Ref**   
**C-LIMS ID**  **Copy Previous**  
**Sample status:** ■

Buildup None   
  Deposition   
  Irradiation  
**Sample Collection Start Date**   
**Sample Collection Stop Date**   
**Sample Ref Date :**   
**Sample Collector :**

mBq   
  Bq   
  kBq  
**Nuclide Vector :**   
**Sample has Certificate**  # additional certificates:

**Vial Type :**   
**Matrix Name :**


| Sample Parameters         | VALUE                               | ERROR                             | UNIT                            |
|---------------------------|-------------------------------------|-----------------------------------|---------------------------------|
| Sample net weight / error | 302.98                              | 0.00                              | (g)                             |
| Sample height:            | 64.28                               | 0.00                              | (mm)                            |
| Additive weight / error   | 0.00                                | 0.00                              | (g)                             |
| Sample Density / error    | 0.54                                | 0.00                              | (g/cm <sup>3</sup> )            |
| Sample quantity / error   | <input type="text" value="0.30"/> C | <input type="text" value="0.00"/> | <input type="text" value="kg"/> |

**Type :**   
**Time Preset (s)**

**Report due date :**  +0 +3 +10 +28  
**Sample disposal date:**

☠ **sample\_risk\_code:**   
**sample\_analysis\_code:**

**500463040020386**

 **Print Sample Barcode**


---


|  | Actual values:                             | Reference values:                          | Difference (%):                        |   |
|--|--|--|--|---|
| <b>Filling height (mm):</b>                  | <input type="text" value="58.68"/>         | <input type="text" value="52.80"/>         | <input type="text" value="10.02"/>     | <b>Requires Efficiency transfert</b> <input checked="" type="checkbox"/>                            |
| <b>Sample Material:</b>                      | <input type="text" value="Dirt1"/>         | <input type="text" value="water"/>         | <input type="text" value="DIFFERENT"/> |   |
| <b>Sample Density: (g/cm<sup>3</sup>)</b>    | <input type="text" value="0.54"/>          | <input type="text" value="1.03"/>          | <input type="text" value="-89.51"/>    | <b>Random Error (%):</b> <input type="text" value="0"/>   |
| <b>Container Diameter: (mm)</b>              | <input type="text" value="113.98"/>        | <input type="text" value="113.98"/>        | <input type="text" value="0.00"/>      | <b>Systematic Error (%):</b> <input type="text" value="0"/>   |
| <b>Container Bottom Thickness: (mm)</b>      | <input type="text" value="1.80"/>          | <input type="text" value="1.80"/>          | <input type="text" value="0.00"/>      | <b>Use fixed Value</b> <input type="checkbox"/>   |
| <b>Container Side Wall: (mm)</b>             | <input type="text" value="2.00"/>          | <input type="text" value="2.00"/>          | <input type="text" value="0.00"/>      | <b>Computed Volume (mL)</b> <input type="text" value="557.45"/>                                     |
| <b>Container Material:</b>                   | <input type="text" value="Polypropyleen"/> | <input type="text" value="Polypropyleen"/> | <input type="text" value="EQUAL"/>     | <b>Reference Volume (mL)</b> <input type="text" value="500.00"/>                                    |
| <b>Container Density: (g/cm<sup>3</sup>)</b> | <input type="text" value="0.91"/>          | <input type="text" value="0.91"/>          | <input type="text" value="0.00"/>      | <b>Filling height / error</b> <input type="text" value="111.14"/> <input type="text" value="0.00"/> |
| <b>Container Gap: (mm)</b>                   | <input type="text" value="3.80"/>          | <input type="text" value="3.80"/>          | <input type="text" value="0.00"/>      |   |




# Efficiency transfer for routine gamma-ray analyses

$$\varepsilon_S(E) = CF(E)\varepsilon_{\text{ref}}(E)$$



- Efficiency transfer  systematic corrections for various counting parameters deviating from reference calibration (EFFTRAN)
- Systematically in each analysis
  - Filling height (if no 100% filling is used)
  - Apparent sample density
  - Sample material
    - **chosen from list** for well known materials
    - Requires specific element composition when gamma-ray energy **(30-100 keV)**
  - Any other sample parameter that is different from the reference


# Efficiency Transfer and Uncertainty Budget imported in Genie 2000



# Uncertainty polynomial


## Known material versus material class

- Known
  - Water
  - PE
  - metals, ...
- Class
  - Organic matter
  - Soil, dirt



- Error curve for each counting geometry
- Filling height, diameter, bodem thickness, sample density, positioning, curve fitting, calibration source, system drift, gamma intensity, counting loses, ...








# Means to determine sample attenuation



- Generic composition
- Comparison with representative standards
  - Not very practical
- XRF-handheld monitor, ICP-MS, ...
  - Specify sample material for the EFFTRAN model >> elements from Mg → U >> ? O, H
- Multi Energy gamma-ray transmission
  - Absolute measurement using collimated beam
  - Relative attenuation measurement using no collimation




# Previous transmission setup



# Calibration curves with transmission measurements



## Ba-133 and Am-241

- X- and  $\gamma$ -rays in the range (21.2 - 303 keV)
- on a set of well known sample materials for which the mass attenuation coefficients are derived from XCOM
- Vial: thickness 10 mm




# Materials


- *Set of materials with well-known chemical composition for which the attenuation data at the transmission energies are computed with XCOM.*
- *Salts, plastics, metals, ...*







# Calibration curve >> E




Transmission source



$$\frac{\mu_{Mat}(E_1)}{\rho_{Mat}} \cdot \frac{\rho_{ref}}{\mu_{ref}(E_1)}$$



$$\frac{\mu_{Mat}(E_n)}{\rho_{Mat}} \cdot \frac{\rho_{ref}}{\mu_{ref}(E_n)}$$





# Calibration curve >> mass att. file

- calibration curves


$$\mu_E \left( \ln \left( \frac{r}{r_{ref}} \right) \right)_{rel} = a_E \ln \left( \frac{r}{r_{ref}} \right) + b_E$$

- $\frac{r}{r_{ref}}(E) \gg$  transmission trough a material relative to the count rate for transmission trough the reference material (water in our case)

- Unknown material X >> mass attenuation coefficients can then be determined from the calibration curves as


$$\mu_X = \mu_0 \frac{\rho_{ref}}{\rho_X} \left( a_E(E) \ln \left( \frac{r_X}{r_{ref}} \right) + b_E \right)$$

→ a mass attenuation file can be set up for an unknown material







# Efficiency Transfer procedure with measured $\mu$ -data




# Transmission spectrum



- 20 keV X-ray complex of Am-241 through a PTFE sample ( $\rho = 2.2 \text{ g/cm}^2$ )
- Attention for a good fit!



- K-edges, a source of discontinuities in gamma attenuation



*Calibration curve for determining  $\mu_{rel}$  at 21.6 keV from transmission relative to water. The materials air, water and aluminium are indicated.  $\mu/\mu_{rel} = 1$  means  $\mu$  is equal to the mass attenuation of water.*


# Implementation routine

2 transmission measurements


Peak surface

Transfer to mass attenuation file > .txt

| TRANSMISSION CALCULATIONS |                     |                           |                                 |  |
|---------------------------|---------------------|---------------------------|---------------------------------|--|
| <b>Parameters</b>         |                     |                           |                                 |  |
| Reference                 | Na2CO3              |                           |                                 |  |
| Date                      | 20/05/2022          |                           |                                 |  |
| Source                    | Am + Ba             |                           |                                 |  |
| Density (g/cm3)           | 1.202               |                           |                                 |  |
| Live time Am241 (s)       | 4201.9              |                           |                                 |  |
| Live time Ba133 (s)       | 4201.9              |                           |                                 |  |
| <b>Source</b>             | <b>Energy (keV)</b> | <b>Peak surface (cts)</b> | <b><math>\mu</math> (cm2/g)</b> | <b>Relative <math>\mu</math> uncertainty</b> |
| Am241                     | 21.16               | 9.30E+04                  | 1.120                           | 0.028  |
| Am241                     | 26.34               | 1.10E+05                  | 0.678                           | 0.272  |
| Ba133                     | 30.97               | 1.30E+06                  | 0.469                           | 0.919  |
| Am241                     | 33.20               | 8.10E+03                  | 0.389                           | 0.258  |
| Ba133                     | 35.05               | 3.54E+05                  | 0.356                           | 1.045  |
| Ba133                     | 53.16               | 3.68E+04                  | 0.161                           | 0.484  |
| Am241                     | 59.54               | 3.21E+06                  | 0.203                           | 0.064  |
| Ba133                     | 81.00               | 6.89E+05                  | 0.133                           | 0.178  |
| Ba133                     | 303.00              | 1.77E+05                  | 0.069                           | 1.180  |



# Conclusion



- The transmission source should be at a distance of at least 20 cm from detector to avoid true coincidence summing effects and measurements are best made with Ba-131 and Am-241 sources sequentially;
- Attenuation parameters can be determined with a precision of a few percent depending on the quality of the calibration curve and the complexity of the spectra;
- No K or L edges in the attenuation of the sample material can be determined with this method;
- The method may fail if spectra are dominated by X-ray fluorescence peaks of the sample material.

# Thank you for your attention

---





## Copyright © SCK CEN

### PLEASE NOTE!

This presentation contains data, information and formats for dedicated use only and may not be communicated, copied, reproduced, distributed or cited without the explicit written permission of SCK CEN.

If this explicit written permission has been obtained, please reference the author, followed by 'by courtesy of SCK CEN'.

Any infringement to this rule is illegal and entitles to claim damages from the infringer, without prejudice to any other right in case of granting a patent or registration in the field of intellectual property.

### **SCK CEN**

Belgian Nuclear Research Centre

Foundation of Public Utility

Registered Office: Avenue Herrmann-Debrouxlaan 40 – BE-1160 BRUSSELS

Operational Office: Boeretang 200 – BE-2400 MOL

# Efftran Transfer data model >> correction function

**EFFTRAN Computation**  
Actual Sample ↔ Reference Sample

E1 → CF(E1)  
E2 → CF(E2)  
E3 → CF(E3)  
E4 → CF(E4)  
....  
En → CF(En)


$$\ln(CF(\ln(E))) = \sum_{k=0}^n C_k \ln(E)^k$$

Spectrum (cam file data)

$$\epsilon_{\text{sam}} = \epsilon_{\text{ref}} \mathbf{CF}$$

$$\ln \epsilon_{\text{ref}}(E) = \sum_{k=0}^n A_k \ln(E)^k$$

$$\ln \epsilon_{\text{corr}}(E) = \sum_{k=0}^n (A_k + C_k) \ln(E)^k$$



Applied to low and high-energy dual efficiency curve


# Cutshall Attenuation Correction Procedures (two times wrong = right ?)

- N.H. Cutshall, correction for  $^{210}\text{Pb}$  in sediments (1981)
- **Near field (wide angle)**


Needs parallel beam of gamma-rays !

Correction only valid for far field geometry !

$$CF(AT) = \frac{-\ln(T)}{1 - T}$$



Bias cutshall vs. EFFTRAN (%) for different geometries



Sediment sample in pillbox (20mm high at different diameters; 50mm, 60mm, 90mm) Cutshall relative to air On a 60mm x 60mm detector But bias is limited!